

---

# **McsPyDataTools Documentation**

***Release 0.3.0***

**Multi Channel Systems MCS GmbH**

**Jun 12, 2018**



CONTENTS

1 McsPyDataTools API Reference 3

1.1 McsPy 3

1.2 The "McsData" module 4

1.3 The "McsCMOS" module 11

2 MCS HDF5 Format Definitions 13

2.1 Definition of the HDF5 format for raw data 13

3 Indices and tables 29

Python Module Index 31



The aim of the McsPyDataTools package is to provide a convenient python interface to access the content of HDF5 data files created by the Multi Channel DataManager (Multi Channel Systems MCS GmbH).

Contents:



## MCSPYDATATOOLS API REFERENCE

### 1.1 McsPy

McsPy is a Python module/package to read, handle and operate on HDF5-based raw data files converted from recordings of devices of the Multi Channel Systems MCS GmbH.

**copyright**

3. 2018 by Multi Channel Systems MCS GmbH

**license** see LICENSE for more details

**class** `McsPy.McsHdf5Protocols`

Class of supported MCS-HDF5 protocol types and version ranges

Entry: (Protocol Type Name => Tuple of supported version range from (including) the first version entry up to (including) the second version entry)

**classmethod** `check_protocol_type_version` (*protocol\_type\_name*, *version*)

Check if the given version of a protocol is supported by the implementation

**Parameters**

- **protocol\_type\_name** – name of the protocol that is tested
- **version** – version number that should be checked

**Returns** is true if the given protocol and version is supported

**class** `McsPy.McsHdf5Types`

Class of supported MCS-HDF5 file structure types and version ranges

Entry: (Protocol TypeID => Tuple of supported version range from (including) the first version entry up to (including) the second version entry)

**classmethod** `check_type_version` (*typeID*, *version*)

Check if the given version of a type is supported by the implementation

**Parameters**

- **protocol\_type\_name** – name of the protocol that is tested
- **version** – version number that should be checked

**Returns** is true if the given type and version is supported

**classmethod** `get_mcs_class_name` (*typeID*, *version=None*)

Returns the McsPy class name, that corresponds to a given Mcs HDF5 file structure type. The function also checks if the requested class supports the Mcs HDF5 file structure type version

**Parameters**

- **typeID** – name of the type that is tested
- **version** – version number that should be checked

**Returns** a McsCMOSMEA class if the given type and version is supported

## 1.2 The “McsData” module

**class** McsPy.McsData.RawData(*raw\_data\_path*)

This class holds the information of a complete MCS raw data file

**recordings**

Access recordings

**class** McsPy.McsData.Recording(*recording\_grp*)

Container class for one recording

Provides the content of the HDF5 *Folder* “Recording\_x” in Python.

**analog\_streams**

Access all analog streams - collection of *AnalogStream* objects

**frame\_streams**

Access all frame streams - collection of *FrameStream* objects

**event\_streams**

Access event streams - collection of *EventStream* objects

**segment\_streams**

Access segment streams - collection of *SegmentStream* objects

**timestamp\_streams**

Access timestamp streams - collection of *TimestampStream* objects

**duration\_time**

Duration of the recording

### 1.2.1 Data-Stream-Structures containing the data

**class** McsPy.McsData.Stream(*stream\_grp*, *info\_type\_name=None*)

Base class for all stream types

**class** McsPy.McsData.AnalogStream(*stream\_grp*)

Container class for one analog stream of several channels. Description for each channel is provided by a channel-associated object of *ChannelInfo*

Provides the content of the HDF5 *Sub-folder* “Stream\_x” of “AnalogStream” in Python.

**get\_channel\_in\_range**(*channel\_id*, *idx\_start*, *idx\_end*)

Get the signal of the given channel over the course of time and in its measured range.

**Parameters**

- **channel\_id** – ID of the channel
- **idx\_start** – index of the first sampled signal value that should be returned (0 <= idx\_start < idx\_end <= count samples)
- **idx\_end** – index of the last sampled signal value that should be returned (0 <= idx\_start < idx\_end <= count samples)



**Returns** Tuple (vector of the signal, unit of the values)

**get\_channel\_sample\_timestamps** (*channel\_id, idx\_start, idx\_end*)

Get the timestamps of the sampled values.

**Parameters**

- **channel\_id** – ID of the channel
- **idx\_start** – index of the first signal timestamp that should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count samples}$ )
- **idx\_end** – index of the last signal timestamp that should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count samples}$ )

**Returns** Tuple (vector of the timestamps, unit of the timestamps)

**class** McsPy.McsData.**FrameStream** (*stream\_grp*)

Container class for one frame stream with different entities

Provides the content of the HDF5 *Subfolder* “*Stream\_x*” of “*FrameStream*” in Python.

**class** McsPy.McsData.**FrameEntity** (*frame\_entity\_group, frame\_info*)

Contains the stream of a specific frame entity. Meta-Information for this entity is available via an associated object of `FrameEntityInfo`

Provides the content of the HDF5 *Subfolder* “*Stream\_x*” of “*FrameStream*” and *Subfolder* “*FrameDataEntity\_x*” in Python.

**get\_sensor\_signal** (*sensor\_x, sensor\_y, idx\_start, idx\_end*)

Get the signal of a single sensor over the course of time and in its measured range.

**Parameters**

- **sensor\_x** – x coordinate of the sensor
- **sensor\_y** – y coordinate of the sensor
- **idx\_start** – index of the first sampled frame that should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count frames}$ )
- **idx\_end** – index of the last sampled frame that should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count frames}$ )

**Returns** Tuple (vector of the signal, unit of the values)

**get\_frame\_timestamps** (*idx\_start, idx\_end*)

Get the timestamps of the sampled frames.

**Parameters**

- **idx\_start** – index of the first sampled frame that should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count frames}$ )
- **idx\_end** – index of the last sampled frame that should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count frames}$ )

**Returns** Tuple (vector of the timestamps, unit of the timestamps)

**class** McsPy.McsData.**EventStream** (*stream\_grp*)

Container class for one event stream with different entities

Provides the content of the HDF5 *Subfolder* “*Stream\_x*” of “*EventStream*” in Python.

**class** McsPy.McsData.**EventEntity** (*event\_data, event\_info*)

Contains the event data of a specific entity. Meta-Information for this entity is available via an associated object of *EventEntityInfo*

Maps data event entity content of the HDF5 *Subfolder* “*Stream\_x*” of “*EventStream*” to Python structures.

**count**

Number of contained events

**get\_events** (*idx\_start=None, idx\_end=None*)

Get all n events of this entity of the given index range ( $\text{idx\_start} \leq \text{idx} < \text{idx\_end}$ )

**Parameters**

- **idx\_start** – start index of the range (including), if nothing is given -> 0
- **idx\_end** – end index of the range (excluding), if nothing is given -> last index

**Returns** Tuple of (2 x n matrix of timestamp (1. row) and duration (2. row), Used unit of time)

**get\_event\_timestamps** (*idx\_start=None, idx\_end=None*)

Get all n event timestamps of this entity of the given index range

**Parameters**

- **idx\_start** – start index of the range, if nothing is given -> 0
- **idx\_end** – end index of the range, if nothing is given -> last index

**Returns** Tuple of (n-length array of timestamps, Used unit of time)

**get\_event\_durations** (*idx\_start=None, idx\_end=None*)

Get all n event durations of this entity of the given index range

**Parameters**

- **idx\_start** – start index of the range, if nothing is given -> 0
- **idx\_end** – end index of the range, if nothing is given -> last index

**Returns** Tuple of (n-length array of duration, Used unit of time)

**class** McsPy.McsData.**SegmentStream** (*stream\_grp*)

Container class for one segment stream of different segment entities

Provides the content of the HDF5 *Subfolder* “*Stream\_x*” of “*SegmentStream*” in Python.

**class** McsPy.McsData.**SegmentEntity** (*segment\_data, segment\_ts, segment\_info*)

Segment entity class, Meta-Information for this entity is available via an associated object of *SegmentEntityInfo*

**DataSybType != Average** → Maps segment entity content of the HDF5 *Subfolder* “*Stream\_x*” of “*SegmentStream*” to Python structures.

**segment\_sample\_count**

Number of contained samples of segments (2d) or multi-segments (3d)

**segment\_count**

Number of segments that are sampled for one time point (2d) -> 1 and (3d) -> n

**get\_segment\_in\_range** (*segment\_id, flat=False, idx\_start=None, idx\_end=None*)

Get the a/the segment signals in its measured range.

**Parameters**

- **segment\_id** – id resp. number of the segment (0 if only one segment is present or the index inside the multi-segment collection)

- **flat** – true -> one-dimensional vector of the sequentially ordered segments, false -> k x n matrix of the n segments of k sample points
- **idx\_start** – index of the first segment that should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count segments}$ )
- **idx\_end** – index of the last segment that should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count segments}$ )

**Returns** Tuple (of a flat vector of the sequentially ordered segments or a k x n matrix of the n segments of k sample points depending on the value of *flat* , and the unit of the values)

**get\_segment\_sample\_timestamps** (*segment\_id*, *flat=False*, *idx\_start=None*, *idx\_end=None*)  
Get the timestamps of the sample points of the measured segment.

#### Parameters

- **segment\_id** – id resp. number of the segment (0 if only one segment is present or the index inside the multi-segment collection)
- **flat** – true -> one-dimensional vector of the sequentially ordered segment timestamps, false -> k x n matrix of the k timestamps of n segments
- **idx\_start** – index of the first segment for that timestamps should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count segments}$ )
- **idx\_end** – index of the last segment for that timestamps should be returned ( $0 \leq \text{idx\_start} < \text{idx\_end} \leq \text{count segments}$ )

**Returns** Tuple (of a flat vector of the sequentially ordered segments or a k x n matrix of the n segments of k sample points depending on the value of *flat* , and the unit of the values)

**class** McsPy.McsData.AverageSegmentTuple (*mean*, *std\_dev*, *time\_tick\_unit*, *signal\_unit*)  
Named tuple that describe one or more average segments (mean, std\_dev, time\_tick\_unit, signal\_unit).

---

#### Note:

- *mean* - mean signal values
  - *std\_dev* - standard deviation of the signal value (it is 0 if there was only one sample segment)
  - *time\_tick\_unit* - sampling interval with time unit
  - *signal\_unit* - measured unit of the signal
- 

#### mean

Alias for field number 0

#### signal\_unit

Alias for field number 3

#### std\_dev

Alias for field number 1

#### time\_tick\_unit

Alias for field number 2

**class** McsPy.McsData.AverageSegmentEntity (*segment\_average\_data*, *segment\_average\_annotation*, *segment\_info*, *seg-*

Contains a number of signal segments that are calculated as averages of number of segments occurred in a given time range. Meta-Information for this entity is available via an associated object of *SegmentEntityInfo*

**DataSybType == Average** → Maps segment entity content of the HDF5 *DataSubType-Average: Subfolder "Stream\_x" of "SegmentStream"* to Python structures.

**number\_of\_averages**

Number of average segments inside this average entity

**sample\_length**

Number of sample points of an average segment

**time\_ranges()**

List of time range tuples for all contained average segments

**Returns** List of tuple with start and end time point

**time\_range(average\_segment\_idx)**

Get the time range for that the average segment was calculated

**Parameters** **average\_segment\_idx** – index resp. number of the average segment

**Returns** Tuple with start and end time point

**average\_counts()**

List of counts of samples for all contained average segments

**Parameters** **average\_segment\_idx** – id resp. number of the average segment

**Returns** sample count

**average\_count(average\_segment\_idx)**

Count of samples that were used to calculate the average

**Parameters** **average\_segment\_idx** – id resp. number of the average segment

**Returns** sample count

**get\_scaled\_average\_segments()**

Get all contained average segments in its measured physical range.

**Returns** *AverageSegmentTuple* containing the k x n matrices for mean and standard deviation of all contained average segments n with the associated sampling and measuring information

**get\_scaled\_average\_segment(average\_segment\_idx)**

Get the selected average segment in its measured physical range.

**Parameters** **segment\_idx** – index resp. number of the average segment

**Returns** *AverageSegmentTuple* containing the mean and standard deviation vector of the average segment with the associated sampling and measuring information

**get\_average\_segments()**

Get all contained average segments AD-offset in ADC values with its measuring conditions

**Returns** *AverageSegmentTuple* containing the mean and standard deviation vector of the average segment in ADC steps with sampling tick and ADC-Step definition

**get\_average\_segment(average\_segment\_idx)**

Get the AD-offset corrected average segment in ADC values with its measuring conditions

**Parameters** **segment\_id** – id resp. number of the segment

**Returns** *AverageSegmentTuple* containing the k x n matrices for mean and standard deviation of all contained average segments in ADC steps with sampling tick and ADC-Step definition

**class** McsPy.McsData.TimestampStream(*stream\_grp*)

Container class for one timestamp stream with different entities

Provides the content of the HDF5 *Subfolder* “Stream\_x” of “TimestampStream” in Python.

**class** McsPy.McsData.TimestampEntity(*timestamp\_data, timestamp\_info*)

Time-Stamp entity class, Meta-Information for this entity is available via an associated object of TimestampEntityInfo

Maps data timestamp entity data of the HDF5 *Subfolder* “Stream\_x” of “TimestampStream” to Python structures.

**count**

Number of contained timestamps

**get\_timestamps** (*idx\_start=None, idx\_end=None*)

Get all n time stamps of this entity of the given index range (*idx\_start* <= *idx* < *idx\_end*)

**Parameters**

- **idx\_start** – start index of the range (including), if nothing is given -> 0
- **idx\_end** – end index of the range (excluding, if nothing is given -> last index

**Returns** Tuple of (n-length array of timestamps, Used unit of time)

## 1.2.2 Info-Classes containing Meta-Information for the data

**class** McsPy.McsData.Info(*info\_data*)

Base class of all info classes

Derived classes contain meta information for data structures and fields.

**get\_field** (*name*)

Get the field with that name -> access to the raw info array

**group\_id**

Get the id of the group that the objects belongs to

**label**

Label of this object

**data\_type**

Raw data type of this object

**class** McsPy.McsData.ChannelInfo(*info\_version, info*)

Contains all describing meta data for one sampled channel

**channel\_id**

Get the ID of the channel

**row\_index**

Get the index of the row that contains the associated channel data inside the data matrix

**adc\_step**

Size and unit of one ADC step for this channel

**version**

Version number of the Type-Definition

**class** McsPy.McsData.InfoSampledData(*info*)

Base class of all info classes for evenly sampled data

**sampling\_frequency**

Get the used sampling frequency in Hz

**sampling\_tick**

Get the used sampling tick

**class** McsPy.McsData.**EventEntityInfo** (*info\_version, info*)

Contains all meta data for one event entity

**id**

Event ID

**raw\_data\_bytes**

Lenght of raw data in bytes

**source\_channel\_ids**

ID's of all channels that were involved in the event generation.

**source\_channel\_labels**

Labels of the channels that were involved in the event generation.

**version**

Version number of the Type-Definition

**class** McsPy.McsData.**SegmentEntityInfo** (*info\_version, info, source\_channel\_infos*)

Contains all meta data for one segment entity

**id**

Segment ID

**pre\_interval**

Interval [start of the segment <- defining event timestamp]

**post\_interval**

Interval [defining event timestamp -> end of the segment]

**type**

Type of the segment like 'Average' or 'Cutout'

**count**

Count of segments inside the segment entity

**version**

Version number of the Type-Definition

**class** McsPy.McsData.**TimeStampEntityInfo** (*info\_version, info*)

Contains all meta data for one timestamp entity

**id**

Timestamp entity ID

**unit**

Unit in which the timestamps are measured

**exponent**

Exponent for the unit in which the timestamps are measured

**measuring\_unit**

Unit in which the timestamp entity was measured

**data\_type**

DataType for the timestamps

**source\_channel\_ids**

ID's of all channels that were involved in the timestamp generation.

**source\_channel\_labels**

Labels of the channels that were involved in the timestamp generation.

**version**

Version number of the Type-Definition

## 1.3 The “McsCMOS” module

### 1.3.1 McsCMOS

Wrapper and Helper to access MCS CMOS Data within H5 Files

**copyright**

3. 2018 by Multi Channel Systems MCS GmbH

**license** see LICENSE for more details

**class** McsPy.McsCMOS.CMOSConvProxy (*parent*)

Private Class, should be embedded within a CMOSData Object. A proxy that transparently converts raw data to calibrated data.

**shape**

Shape of the data

**class** McsPy.McsCMOS.CMOSData (*path*)

Wrapper for a HDF5 File containing CMOS Data

**class** McsPy.McsCMOS.CMOSSpikes (*path*)

Wrapper for a HDF5 File containing CMOS Spike Data. Spike Information is accessible through the .spike Member, Waveform Information (if available) through the .waveforms Member.





## MCS HDF5 FORMAT DEFINITIONS

### 2.1 Definition of the HDF5 format for raw data

MCS-HDF5 Protocol Type: RawData (Raw-Data protocol)

Protocol Version: 3 based on the definitions of RawDataFileIO in version 10.

All strings are only ASCII-encoded

#### 2.1.1 Changelog

Version 1:

- Initial draft

Version 2:

- New Root-Folder attributes added to detect name and version of the creating application and library

Version 3:


- Data structures for **DataSubType::Average** of **StreamType::Segment** added

#### Hierarchy

##### Root-Folder “/”

Contains all information for one experiment - measured data (inside the folder Data) and a description (possibly in the future) inside the folder Experiment/Description/...

Attributes:



Name	Description	Data Type	MCS-HDF5 Protocol Version
McsHdf5ProtocolType	Type of the used MCS-HDF5 protocol definition (e.g. <b>RawData</b> for the raw data MCS-HDF5 definitions)	[String,Scalar]	1 ≤
McsHdf5ProtocolVersion	Version number of the used MCS-HDF5 protocol	[Integer,Scalar]	1 ≤
GeneratingApplicationName	Name of the application that generated this HDF5 file	[String,Scalar]	2 ≤
GeneratingApplicationVersion	Version of the application that generated this HDF5 file	[String,Scalar]	2 ≤
McsDataToolsVersion	Version of the McsDataTools library that was used by the application to create the HDF5 file	[String,Scalar]	2 ≤

Datasets:

- none

### Folder “Data”

Navigation: /Data

Contains all recordings for this experiment.

Attributes:

Name	Description	Data Type
ProgramName	Name of the recording program	[String,Scalar]
ProgramVersion	Version number of the recording program	[String,Scalar]
MeaName	Name of the recorded MEA	[String,Scalar]
MeaLayout	Layout descriptor	[String,Scalar]
MeaSN	Serial number of the MEA	[String,Scalar]
Date	Date of the recording	[String,Scalar]
DateInTicks	Date of the recording in .NET ticks (100 ns)	[Long(64-bit Integer),Scalar]
FileGUID	GUID of the converted raw data file	[String,Scalar]
Comment	Comment	[String,Scalar]

Datasets:

- none

### Folder “Recording\_x”

Navigation: /Data/Recording\_x

Contains all recorded streams for recording x.

Attributes:

Name	Description	Data Type
RecordingID	Recording ID	[Integer(32-bit Integer),Scalar]
RecordingType	Recording type	[String,Scalar]
TimeStamp	Start time of the recording in microseconds	[Long(64-bit Integer),Scalar]
Duration	Total recording duration in microseconds (This duration can differ from the actual duration of the recorded data!!!)	[Long(64-bit Integer),Scalar]
Label	Label	[String,Scalar]
Comment	Comment	[String,Scalar]

Datasets:

- none

### Folder “AnalogStream”

Navigation: /Data/Recording\_x/AnalogStream

(Organisational) folder for all channel-based streams of this recording

Attributes:

- none

Datasets:

- none

### Sub-folder “Stream\_x” of “AnalogStream”

Navigation: /Data/Recording\_x/AnalogStream/Stream\_x

Container for an analog stream

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[ <b>Int</b> (32-bit Integer), <b>Scalar</b> ]	$1 \leq$
Label	Label	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
SourceStreamGUID	GUID of the source streams	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
StreamGUID	GUID	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
StreamType	Type of the stream, e.g. <b>Electrode</b>	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
DataSubType	Sub-type of the analog stream (e.g. <b>Analog</b> )	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$

Datasets:

- Matrix InfoChannel  $\rightarrow n \times 16$  matrix of describing information vectors for the  $n$  channels:
  - Attributes: InfoVersion  $\rightarrow$  Version number of the Info-Objects [**Int**(32-bit Integer),**Scalar**]

Name	Description	Data Type	InfoVersion
ChannelID	ID of the channel as given by the recording software	[Int(32-bit Integer),Array(Size 1)]	1 ≤
RowIndex	Row number of this channel inside the <b>ChannelData</b> matrix where the data of this channel is stored	[Int(32-bit Integer),Array(Size 1)]	1 ≤
GroupID	ID of the group that this channel belongs to	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Label	Label of the channel	[String,Array]	1 ≤
RawDataType	Type of the raw data	[String,Array]	1 ≤
Unit	Physical unit of the measured sensor value	[String,Array]	1 ≤
Exponent	Exponent $n \Rightarrow 10^n$ in which the channel values magnitude is measured (e.g. k,m, $\mu$ ,...)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
ADZero	ADC-Step that represents the 0-point of the measuring range of the ADC	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Tick	Sample tick $\Delta$ between two sample points of a channel in $\mu$ s $\Rightarrow$ sampling frequency = $1000000 / \Delta$	[Long(64-bit Integer),Array(Size 1)]	1 ≤
ConversionFactor	Conversion factor for the mapping ADC-Step $\Rightarrow$ measured value	[Long(64-bit Integer),Array(Size 1)]	1 ≤
ADCBits	Number of bits used by the AD-Converter	[Int(32-bit Integer),Array(Size 1)]	1 ≤
HighPassFilterType	Type of the high-pass filter (empty string if not available)	[String,Scalar]	1 ≤
HighPassFilterCutOffFrequency	Cut-off frequency of the high-pass filter ('-1'-String if not available)	[String,Scalar]	1 ≤
HighPassFilterOrder	Order of the high-pass filter (-1 if not available)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
LowPassFilterType	Type of the low-pass filter (empty string if not available)	[String,Scalar]	1 ≤
LowPassFilterCutOffFrequency	Cut-off frequency of the low-pass filter ('-1'-String if not available)	[String,Scalar]	1 ≤
LowPassFilterOrder	Order of the low-pass filter (-1 if not available)	[Int(32-bit Integer),Array(Size 1)]	1 ≤

- 2-dimensional Data-Matrix `ChannelData` → Data for sampled channels organized as  $n \times m$  matrix  $\Rightarrow$  one row per channel and one column per sample time point
  - reconstruct the value of the measured signal:  $y(\text{channel}, t_{ind}) = (\text{ChannelData}[\text{InfoChannel}[\text{channel}].\text{RowIndex}, t_{ind}] - \text{ADZero}) * \text{InfoChannel}[\text{channel}].\text{ConversionFactor} * 10^{\text{InfoChannel}[\text{channel}].\text{Exponent}}$  in `InfoChannel[channel].Unit`
  - reconstruct the sample time point:  $t = t_{ind} * \text{InfoChannel}[\text{channel}].\text{Tick}$  in  $\mu s$
- Matrix `ChannelDataTimeStamps` →  $k \times 3$  matrix of segments where the rows are one segment and the columns are:
  - first column → time stamp of the first sample point of the segment
  - second column → first index (column) of the segment in **ChannelData**
  - third column → last index (column) of the segment in **ChannelData**

### Folder “FrameStream”

Navigation: `/Data/Recording_x/FrameStream`

(Organisational) folder for all frame-based streams of this recording

Attributes:

- none

Datasets:

- none

### Subfolder “Stream\_x” of “FrameStream”

Navigation: `/Data/Recording_x/FrameStream/Stream_x`

Folder that contains all Frame-Entities of one Frame-Stream:

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[ <b>Int</b> (32-bit Integer), <b>Scalar</b> ]	$1 \leq$
Label	Label	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
SourceStreamGUID	GUID of the source stream	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
StreamGUID	GUID	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
StreamType	Type of the stream <b>Frame</b>	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
DataSubType	Sub-type of the event stream (e.g. <b>Spike-TimeStamp</b> )	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$

Datasets:

- Matrix `InfoFrame` →  $n \times 24$  matrix of describing information vectors for the  $n$  Frame-Entities:
  - Attributes: `InfoVersion` → Version number of the Info-Objects [**Int**(32-bit Integer), **Scalar**]

Name	Description	Data Type	InfoVersion
FrameID	ID of the frame entity as given by the recording software	[Int(32-bit Integer),Array(Size 1)]	1 ≤
FrameDataID	ID of the frame entity inside the stream folder that maps this information vector to the entity folder ( <b>FrameDataID</b> → subfolder <b>FrameDataEntity_FrameDataID</b> )	[Int(32-bit Integer),Array(Size 1)]	1 ≤
GroupID	ID of the group that this frame entity belongs to	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Label	Label of the entity	[String,Array]	1 ≤
RawDataType	Type of the raw data	[String,Array]	1 ≤
Unit	Physical unit of the measured sensor value	[String,Array]	1 ≤
Exponent	Exponent $n \Rightarrow 10^n$ in which the sensor values magnitude is measured (e.g. k,m, $\mu$ ,...)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
ADZero	ADC-Step that represents the 0-point of the measuring range of the ADC	[Int(32-bit Integer),Array(Size 1)]	1 ≤
ADCBits	Number of bits used by the AD-Converter	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Tick	Sample tick $\Delta$ between two frames in $\mu$ s $\Rightarrow$ sampling frequency = $1000000 / \Delta$	[Long(64-bit Integer),Array(Size 1)]	1 ≤
HighPassFilterType	Type of the high-pass filter (empty string if not available)	[String,Scalar]	1 ≤
HighPassFilterCutOffFrequency	Cut-off frequency of the high-pass filter ('-1'-String if not available)	[String,Scalar]	1 ≤
HighPassFilterOrder	Order of the high-pass filter (-1 if not available)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
LowPassFilterType	Type of the low-pass filter (empty string if not available)	[String,Scalar]	1 ≤
LowPassFilterCutOffFrequency	Cut-off frequency of the low-pass filter ('-1'-String if not available)	[String,Scalar]	1 ≤
LowPassFilterOrder	Order of the low-pass filter (-1 if not available)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
SensorSpacing	Distance between adjacent sensors in $\mu$ m	[Int(32-bit Integer),Array(Size 1)]	1 ≤



### Subfolder “FrameDataEntity\_x”

Navigation: /Data/Recording\_x/FrameStream/Stream\_x/FrameDataEntity\_x

Contains all datasets of the Frame-Entity x

Datasets:

- Matrix ConversionFactors  $\rightarrow n \times m$  matrix of conversion factors for the sensor array
- 3-dimensional Data-Cube FrameData  $\rightarrow$  cube of the frame data organized as one frame to one sample time point ( $n \times m$  matrix of sampled signal values per sensor)  $\times$  sample time points
  - reconstruct the value of the measured signal:  $y = (\text{FrameData}[x,y,t] - \text{ADZero}) * \text{ConversionFactors}[x,y]$
  - reconstruct the sample time point:
- Matrix FrameDataTimeStamps  $\rightarrow k \times 3$  matrix of segments where the rows are one segment and the columns are:
  - first column  $\rightarrow$  time stamp of the first sample point of the segment
  - second column  $\rightarrow$  first index (z-axis) of the segment in **FrameData**
  - third column  $\rightarrow$  last index (z-axis) of the segment in **FrameData**

Datasets:

- none

### Folder “EventStream”

Navigation: /Data/Recording\_x/EventStream

(Organisational) folder for all event-based streams of this recording

Attributes:

- none

Datasets:

- none

### Subfolder “Stream\_x” of “EventStream”

Navigation: /Data/Recording\_x/EventStream/Stream\_x

Folder that contains all Event-Entities of one Event-Stream:

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[ <b>Int</b> (32-bit Integer), <b>Scalar</b> ]	$1 \leq$
Label	Label	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
SourceStreamGUID	GUID of the source stream	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
StreamGUID	GUID of the current stream	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
StreamType	Type of the stream <b>Event</b>	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$
DataSubType	Sub-type of the event stream (e.g. <b>StgSideband</b> , <b>UserInput</b> , <b>DigitalPort</b> )	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$

Sub-type Description:

- **StgSideband** → The event is associated to a STG sideband change.
- **UserInput** → The event is associated with an user input.
- **DigitalPort** → The event is associated with a digital port change.

Datasets:

- Matrix **InfoEvent** →  $n \times 7$  matrix of describing information vectors for the  $n$  Event-Entities:
  - Attributes: **InfoVersion** → Version number of the Info-Objects [**Int**(32-bit Integer),**Scalar**]

Name	Description	Data Type	InfoVersion
EventID	ID of the event entity	[ <b>Int</b> (32-bit Integer), <b>Array</b> (Size 1)]	$1 \leq$
GroupID	ID of the group that the entity belongs to	[ <b>Int</b> (32-bit Integer), <b>Array</b> (Size 1)]	$1 \leq$
Label	Label of the entity	[ <b>String</b> , <b>Array</b> ]	$1 \leq$
RawDataType	Type of the raw data	[ <b>String</b> , <b>Array</b> ]	$1 \leq$
RawDataBytes	Number of bytes of the raw data type	[ <b>Int</b> (32-bit Integer), <b>Array</b> (Size 1)]	$1 \leq$
SourceChannelIDs	Comma separated list of ID's of (source) channel that were involved in the generation of this event	[ <b>String</b> , <b>Array</b> ]	$1 \leq$
SourceChannelLabels	Comma separated list of labels of the source channels	[ <b>String</b> , <b>Scalar</b> ]	$1 \leq$

- 2-dimensional matrix **EventEntity\_x** →  $2 \times n$  matrix  $\Rightarrow$   $n$  events with describing vector (time stamp of event, duration of event)
  - Attributes: Short description of content
  - $t_{\text{event } i} = \text{EventEntity\_x}[0, i]$  in  $\mu s$
  - $\Delta_{\text{event } i} = \text{EventEntity\_x}[1, i]$  in  $\mu s$

## Folder “SegmentStream”

Navigation: /Data/Recording\_x/SegmentStream

(Organisational) folder for all segment-based streams of this recording. A segment is a cutout of parts of the sampled signal relative to an event, defined by a pre- and post interval.

Attributes:

- none

Datasets:

- none

## Subfolder “Stream\_x” of “SegmentStream”

Navigation: /Data/Recording\_x/SegmentStream/Stream\_x

Folder that contains all Segment-Entities of one Segment-Stream:

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[ <b>Int</b> (32-bit Integer), <b>Scalar</b> ]	1 ≤
Label	Label	[ <b>String</b> , <b>Scalar</b> ]	1 ≤
SourceStreamGUID	GUID of the source stream	[ <b>String</b> , <b>Scalar</b> ]	1 ≤
StreamGUID	GUID of the current stream	[ <b>String</b> , <b>Scalar</b> ]	1 ≤
StreamType	Type of the stream <b>Segment</b>	[ <b>String</b> , <b>Scalar</b> ]	1 ≤
DataSubType	Sub-type of the segment stream (e.g. <b>Spike</b> )	[ <b>String</b> , <b>Scalar</b> ]	1 ≤

Datasets:

- Matrix InfoSegment  $\rightarrow n \times 7$  matrix of describing information vectors for the n Segment-Entities:
  - Attributes: InfoVersion  $\rightarrow$  Version number of the Info-Objects [**Int**(32-bit Integer),**Scalar**]

Name	Description	Data Type	InfoVersion
SegmentID	ID of the segment entity	[ <b>Int</b> (32-bit Integer), <b>Array</b> (Size 1)]	1 ≤
GroupID	ID of the group that the segment entity belongs to	[ <b>Int</b> (32-bit Integer), <b>Array</b> (Size 1)]	1 ≤
Label	Label of the entity	[ <b>String</b> , <b>Array</b> ]	1 ≤
PreInterval	Time interval in $\mu$ s before the segment defining event occurred - definition of the beginning of the segment	[ <b>Int</b> (64-bit Integer), <b>Array</b> (Size 1)]	1 ≤
PostInterval	Time interval in $\mu$ s after the segment defining event occurred - definition of the end of the segment  length of the segment = PreInterval + PostInterval in $\mu$ s	[ <b>Int</b> (64-bit Integer), <b>Array</b> (Size 1)]	1 ≤
SegmentType	Type of the segment (e.g. <b>SpikeCutout</b> )	[ <b>String</b> , <b>Array</b> ]	1 ≤
SourceChannelIDs	Comma separated list of ID's of (source) channels that the segments are taken from → Link to the SourceChannelInfo matrix	[ <b>String</b> , <b>Array</b> ]	1 ≤

- 2-dimensional matrix SourceChannelInfo  $\rightarrow n \times 15$  matrix  $\Rightarrow n$  of describing vectors for the  $n$  source channels, the structure is the same as in ChannelInfo used in section *Sub-folder "Stream\_x" of "AnalogStream"*
  - Attributes: InfoVersion  $\rightarrow$  Version number of the Info-Objects [**Int**(32-bit Integer),**Scalar**]
- Vector SegmentData\_ts\_x  $\rightarrow n$  time stamps in  $\mu$ s of the event triggering the segment, one for each of the  $n$  segments contained by segment entity  $x$
- 2-dimensional matrix or 3-dimensional cube SegmentData\_x  $\rightarrow k \times n$  matrix ( $k$  sample points for one segment,  $n$  number of sampled segments) or  $k \times m \times n$  cube ( $k$  sample points for one segment,  $m$  number of segments for one time point/for one multi-segment,  $n$  number of sampled multi-segments) of segment data:
  - Attributes: SourceChannelID  $\rightarrow$  Comma separated list of ID's of (source) channels that the segments are taken from [**String**,**Scalar**] (the same as in InfoSegment, repeated for clarification)
  - reconstruct the value of the measured segment signal (only one segment  $id_{\text{segment}} \rightarrow$  2-dimensional matrix

$M[\text{row}, \text{col}]$ :

$$* \quad t_{ind}[\text{row}, \text{col}] = \text{SegmentData ts } x[\text{col}] + (\text{row} - 1) * \text{tick}_{\text{source-channel}} - \text{PreInterval in } \mu\text{s}$$

$$* \quad y(id_{\text{segment}}, t_{ind}(\text{row}, \text{col})) = \frac{(\text{SegmentData } x[\text{row}, \text{col}] - \text{ADZero}_{\text{source-channel}})}{\text{ConversionFactor}_{\text{source-channel}} * 10^{\text{Exponent}_{\text{source-channel}}}} \text{ in InfoChannel}[\text{source-channel}].\text{Unit}$$

- reconstruct the value of the measured segment signal (m segments  $\rightarrow$  multi-segments  $\rightarrow$  3-dimensional cube  $M[\text{row}, \text{col}, z]$ ):

$$* \quad \text{col} \rightarrow id_{\text{segment}} \rightarrow \text{source-channel}$$

$$* \quad t_{ind}[\text{row}, \text{col}, z] = \text{SegmentData ts } x[z] + (\text{row} - 1) * \text{tick}_{\text{source-channel}[\text{col}]} \text{ in } \mu\text{s}$$

$$* \quad y(id_{\text{segment}}, t_{ind}(\text{row}, z)) = \frac{(\text{SegmentData } x[\text{row}, \text{col}, z] - \text{ADZero}_{\text{source-channel}[\text{col}]})}{\text{ConversionFactor}_{\text{source-channel}[\text{col}]} * 10^{\text{Exponent}_{\text{source-channel}[\text{col}]}}} \text{ in InfoChannel}[\text{source-channel}[\text{col}]].\text{Unit}$$

### DataSubType-Average: Subfolder “Stream\_x” of “SegmentStream”

Navigation: /Data/Recording\_x/SegmentStream/Stream\_x

Folder that contains all Segment-Entities of one Segment-Stream with **DataSybType == Average**:

Attributes: no difference to the standard case above

Datasets:

- Matrix InfoSegment: no difference to the standard case above
- Matrix SourceChannelInfo: no difference to the standard case above
- $(3 \times n)$  matrix AverageData\_Range\_x  $\rightarrow$  (**start**, **end**, **count**) per segment average  $\times$  count of segment averages contained by segment entity x. **start** and **end** denote the start and end timestamp in  $\mu\text{s}$  of the interval that contains all averaged segments. **count** is the number of averaged segments.
  - Attributes: description of the content
- $(2 \times k \times n)$  cube AverageData\_x  $\rightarrow$  (mean and standard deviation)  $\times$  k sample points of the segment  $\times$  n number of segment averages
  - Attributes:: description of the content
  - reconstruct the value of the mean and standard deviation of the average segment (n average segments  $\rightarrow$  3-dimensional cube  $M[\text{row}, \text{col}, z]$ ):
    - \* row: mean  $\rightarrow$  row = 0; StdDev  $\rightarrow$  row = 1
    - \* col:  $t_{ind}(\text{col}) = (\text{col} - 1) * \text{tick}_{\text{source-channel}} \rightarrow$  time range  $(0, \text{PreInterval}[\text{SegmentID}] + \text{PreInterval}[\text{SegmentID}])$  in  $\mu\text{s}$
    - \* z:  $z = id_{\text{average}}$  (number of average segment)
    - \*  $Mean(id_{\text{average}}, t_{ind}(\text{col})) = \frac{(\text{AverageData } x[0, \text{col}, id_{\text{average}}] - \text{ADZero}_{\text{source-channel}})}{\text{ConversionFactor}_{\text{source-channel}} * 10^{\text{Exponent}_{\text{source-channel}}}} \text{ in InfoChannel}_{\text{source-channel}}.\text{Unit}$
    - \*  $StdDev(id_{\text{average}}, t_{ind}(\text{col})) = \text{AverageData } x[1, \text{col}, id_{\text{average}}] * \text{ConversionFactor}_{\text{source-channel}} * 10^{\text{Exponent}_{\text{source-channel}}} \text{ in InfoChannel}_{\text{source-channel}}.\text{Unit}$

### Folder “TimeStampStream”

Navigation: /Data/Recording\_x/TimeStampStream

(Organisational) folder for all TimeStamp-based streams of this recording

Attributes:

- none

Datasets:

- none

### Subfolder “Stream\_x” of “TimeStampStream”

Navigation: /Data/Recording\_x/TimeStampStream/Stream\_x

Folder that contains all TimeStamp-Entities of one TimeStamp-Stream:

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[ <b>Int</b> (32-bit Integer), <b>Scalar</b> ]	1 ≤
Label	Label	[ <b>String</b> , <b>Scalar</b> ]	1 ≤
SourceStreamGUID	GUID of the source stream	[ <b>String</b> , <b>Scalar</b> ]	1 ≤
StreamGUID	GUID of the current stream	[ <b>String</b> , <b>Scalar</b> ]	1 ≤
StreamType	Type of the stream <b>TimeStamp</b>	[ <b>String</b> , <b>Scalar</b> ]	1 ≤
DataSubType	Sub-type of the TimeStamp stream (e.g. <b>NeuralSpike</b> )	[ <b>String</b> , <b>Scalar</b> ]	1 ≤

Sub-type Description:

- **NeuralSpike** → The entity contains time stamps of neural spikes

Datasets:

- Matrix InfoTimeStamp →  $n \times 7$  matrix of describing information vectors for the  $n$  Event-Entities:
  - Attributes: InfoVersion → Version number of the Info-Objects [**Int**(32-bit Integer),**Scalar**]

Name	Description	Data Type	InfoVersion
TimeStampEntityID	ID of the event entity	[Int(32-bit Integer),Array(Size 1)]	1 ≤
GroupID	ID of the group that the entity belongs to	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Label	Label of the entity	[String,Array]	1 ≤
Unit	Physical unit of the measured sensor value	[String,Array]	1 ≤
Exponent	Exponent $n \Rightarrow 10^n$ in which the channel values magnitude is measured (e.g. k,m, $\mu$ ,...)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
SourceChannelIDs	Comma separated list of ID's of (source) channel that were involved in the generation of this event	[String,Array]	1 ≤
SourceChannelLabels	Comma separated list of labels of the source channels	[String,Scalar]	1 ≤

- Vector TimeStampEntity\_x  $\rightarrow$  n time stamps in  $\mu s$

### Comment

All time-related information except dates (100ns ticks) are given in  $\mu s$  ticks!!

Category:Software





## INDICES AND TABLES

- `genindex`
- `modindex`
- `search`



## PYTHON MODULE INDEX

### m

McsPy, [3](#)

McsPy.McsCMOS, [11](#)

McsPy.McsData, [4](#)

## INDEX

- adc\_step (McsPy.McsData.ChannelInfo attribute), 9
- analog\_streams (McsPy.McsData.Recording attribute), 4
- AnalogStream (class in McsPy.McsData), 4
- average\_count() (McsPy.McsData.AverageSegmentEntity method), 8
- average\_counts() (McsPy.McsData.AverageSegmentEntity method), 8
- AverageSegmentEntity (class in McsPy.McsData), 7
- AverageSegmentTuple (class in McsPy.McsData), 7
- channel\_id (McsPy.McsData.ChannelInfo attribute), 9
- ChannelInfo (class in McsPy.McsData), 9
- check\_protocol\_type\_version() (McsPy.McsHdf5Protocols class method), 3
- check\_type\_version() (McsPy.McsHdf5Types class method), 3
- CMOSConvProxy (class in McsPy.McsCMOS), 11
- CMOSData (class in McsPy.McsCMOS), 11
- CMOSSpikes (class in McsPy.McsCMOS), 11
- count (McsPy.McsData.EventEntity attribute), 6
- count (McsPy.McsData.SegmentEntityInfo attribute), 10
- count (McsPy.McsData.TimeStampEntity attribute), 9
- data\_type (McsPy.McsData.Info attribute), 9
- data\_type (McsPy.McsData.TimeStampEntityInfo attribute), 10
- duration\_time (McsPy.McsData.Recording attribute), 4
- event\_streams (McsPy.McsData.Recording attribute), 4
- EventEntity (class in McsPy.McsData), 5
- EventEntityInfo (class in McsPy.McsData), 10
- EventStream (class in McsPy.McsData), 5
- exponent (McsPy.McsData.TimeStampEntityInfo attribute), 10
- frame\_streams (McsPy.McsData.Recording attribute), 4
- FrameEntity (class in McsPy.McsData), 5
- FrameStream (class in McsPy.McsData), 5
- get\_average\_segment() (McsPy.McsData.AverageSegmentEntity method), 8
- get\_average\_segments() (McsPy.McsData.AverageSegmentEntity method), 8
- get\_channel\_in\_range() (McsPy.McsData.AnalogStream method), 4
- get\_channel\_sample\_timestamps() (McsPy.McsData.AnalogStream method), 5
- get\_event\_durations() (McsPy.McsData.EventEntity method), 6
- get\_event\_timestamps() (McsPy.McsData.EventEntity method), 6
- get\_events() (McsPy.McsData.EventEntity method), 6
- get\_field() (McsPy.McsData.Info method), 9
- get\_frame\_timestamps() (McsPy.McsData.FrameEntity method), 5
- get\_mcs\_class\_name() (McsPy.McsHdf5Types class method), 3
- get\_scaled\_average\_segment() (McsPy.McsData.AverageSegmentEntity method), 8
- get\_scaled\_average\_segments() (McsPy.McsData.AverageSegmentEntity method), 8
- get\_segment\_in\_range() (McsPy.McsData.SegmentEntity method), 6
- get\_segment\_sample\_timestamps() (McsPy.McsData.SegmentEntity method), 7
- get\_sensor\_signal() (McsPy.McsData.FrameEntity method), 5
- get\_timestamps() (McsPy.McsData.TimeStampEntity method), 9
- group\_id (McsPy.McsData.Info attribute), 9
- id (McsPy.McsData.EventEntityInfo attribute), 10
- id (McsPy.McsData.SegmentEntityInfo attribute), 10
- id (McsPy.McsData.TimeStampEntityInfo attribute), 10
- Info (class in McsPy.McsData), 9
- InfoSampledData (class in McsPy.McsData), 9
- label (McsPy.McsData.Info attribute), 9
- McsHdf5Protocols (class in McsPy), 3
- McsHdf5Types (class in McsPy), 3

[McsPy \(module\), 3](#)  
[McsPy.McsCMOS \(module\), 11](#)  
[McsPy.McsData \(module\), 4](#)  
[mean \(McsPy.McsData.AverageSegmentTuple attribute\), 7](#)  
[measuring\\_unit \(McsPy.McsData.TimestampEntityInfo attribute\), 10](#)  
[number\\_of\\_averages \(McsPy.McsData.AverageSegmentEntity attribute\), 8](#)  
[post\\_interval \(McsPy.McsData.SegmentEntityInfo attribute\), 10](#)  
[pre\\_interval \(McsPy.McsData.SegmentEntityInfo attribute\), 10](#)  
[raw\\_data\\_bytes \(McsPy.McsData.EventEntityInfo attribute\), 10](#)  
[RawData \(class in McsPy.McsData\), 4](#)  
[Recording \(class in McsPy.McsData\), 4](#)  
[recordings \(McsPy.McsData.RawData attribute\), 4](#)  
[row\\_index \(McsPy.McsData.ChannelInfo attribute\), 9](#)  
[sample\\_length \(McsPy.McsData.AverageSegmentEntity attribute\), 8](#)  
[sampling\\_frequency \(McsPy.McsData.InfoSampledData attribute\), 9](#)  
[sampling\\_tick \(McsPy.McsData.InfoSampledData attribute\), 10](#)  
[segment\\_count \(McsPy.McsData.SegmentEntity attribute\), 6](#)  
[segment\\_sample\\_count \(McsPy.McsData.SegmentEntity attribute\), 6](#)  
[segment\\_streams \(McsPy.McsData.Recording attribute\), 4](#)  
[SegmentEntity \(class in McsPy.McsData\), 6](#)  
[SegmentEntityInfo \(class in McsPy.McsData\), 10](#)  
[SegmentStream \(class in McsPy.McsData\), 6](#)  
[shape \(McsPy.McsCMOS.CMOSConvProxy attribute\), 11](#)  
[signal\\_unit \(McsPy.McsData.AverageSegmentTuple attribute\), 7](#)  
[source\\_channel\\_ids \(McsPy.McsData.EventEntityInfo attribute\), 10](#)  
[source\\_channel\\_ids \(McsPy.McsData.TimestampEntityInfo attribute\), 10](#)  
[source\\_channel\\_labels \(McsPy.McsData.EventEntityInfo attribute\), 10](#)  
[source\\_channel\\_labels \(McsPy.McsData.TimestampEntityInfo attribute\), 11](#)  
[std\\_dev \(McsPy.McsData.AverageSegmentTuple attribute\), 7](#)  
[Stream \(class in McsPy.McsData\), 4](#)  
[time\\_range\(\) \(McsPy.McsData.AverageSegmentEntity method\), 8](#)  
[time\\_ranges\(\) \(McsPy.McsData.AverageSegmentEntity method\), 8](#)  
[time\\_tick\\_unit \(McsPy.McsData.AverageSegmentTuple attribute\), 7](#)  
[timestamp\\_streams \(McsPy.McsData.Recording attribute\), 4](#)  
[TimestampEntity \(class in McsPy.McsData\), 9](#)  
[TimestampEntityInfo \(class in McsPy.McsData\), 10](#)  
[TimestampStream \(class in McsPy.McsData\), 8](#)  
[type \(McsPy.McsData.SegmentEntityInfo attribute\), 10](#)  
[unit \(McsPy.McsData.TimestampEntityInfo attribute\), 10](#)  
[version \(McsPy.McsData.ChannelInfo attribute\), 9](#)  
[version \(McsPy.McsData.EventEntityInfo attribute\), 10](#)  
[version \(McsPy.McsData.SegmentEntityInfo attribute\), 10](#)  
[version \(McsPy.McsData.TimestampEntityInfo attribute\), 11](#)